

Medical Master

Electric Hospital Bed MEB-253 (3 motors)

Standard hospital care bed of MEB-253 line provide very good level of comfort for user, caregiver and nursing staff.

This full electric adjustable bed has 4 sections in standard and electric operation for all movements.

Both models could add different options to use for different purposes, like weight scale, X-ray cassette, lifting pole, Balkan frame, battery...etc.

Technical Specification:

Size: L2100 x W1000 mm +/- 5%

Adjustable height: 400~700 mm +/- 5%

Back adjustable: 0~80°

Knee adjustable: 0~35°

Loading capacity: 180 kg (253)

Key Feature:

1. Easy to clean with PP detachable platform with coating steel support (HFB-07)
2. Material in reach of patient care are anti-bacterial (power coating steel)
3. Full electric adjustable movement for easy handling
4. Ergonomic hand controller
5. Removable PP head & foot board with rotating bumper
6. Collapsible 3/4 length side rail
 - Single push button for fold down
 - Aluminum pillar with rail
7. 125 mm high quality brake caster
8. Urinal hooks on both side of bed
9. IV pole socket at 4 corners of bed
10. 10 cm normal mattress

MEB-253:

