



YZ Series

Infant Oxygen Hood



Care with love,
Creation of comfortable oxygen inhalation



Endless Love To Guard The New Life Safe

Effectively avoid oxygen toxicity, reduce the respiratory resistance of newborns, and improve brain hypoxia in the process of oxygen supply. Meet different needs of various neonatal respiratory diseases by adjusting oxygen concentration in the hood.

Oxygen Supply Environment - Safer and More Stable

- Uniform and stable oxygen concentration, can be adjusted as needed.
Aerodynamics designed air-guide valve, making oxygen flow in the hood more uniform. Airflow noise can be eliminated to create a quiet oxygen therapy environment (YZ-260/300). Flexible selection of the number of exposed holes, smooth operation without damping, adjustable air-oxygen ratio (YZ-130/160/200).
- The high flow rate gas can accelerate the exhaled carbon dioxide from newborns to be discharged from the hood and avoid repeated inhalation.
- Reduced irritation to the nasal mucosa.



Clinical Nursing - Simpler and More Efficient YZ-260/YZ-300

- Silicone curtain piece is designed between the hood and the baby's neck for transition, thereby ensuring airtightness and preventing the newborns from being scratched.
- The side door can be slid up and down without moving the hood, thereby facilitating treatment and care by medical staff and always maintaining an oxygen-rich environment.
- Two brackets are provided to meet different sizes of oxygen supply hoses.



Streamlined Design - Smaller and More Flexible

- Streamlined appearance design
Integrated injection molding, fully transparent visualization, allowing medical staff to easily observe the treatment of newborns in the hood at a distance.
Made of PC material, colorless and environmentally friendly, durable, compact and portable, easy to disassemble and clean (YZ-260/300).
- Flexible adaptation according to clinical needs
Five optional specifications including YZ-200 (large), YZ-160 (medium), YZ-130 (small), YZ-260 and YZ-300, which can be flexibly selected according to the size of the head and neck of newborns.



YZ-130/YZ-160/YZ-200



YZ-260



YZ-300

Ningbo David Medical Device Co.,Ltd

Add: Bldg No. 5 Ningbo Smart Park, #98 Chuangyuan Rd, Hi-tech Zone, Ningbo, China

Tel: 0086-574-8780 0008 or 8780 0007

P.C.: 315048

E-mail: sales@nbdavid.com

Fax: 0086-574-8780 3714

Web: www.nbdavid.com/en/

Specifications are changed without notice for performance improvement.

Printed date: May 2022 Version No.: 1

