

exacta®



This classic twin-tube sphygmomanometer has been further improved and reduced in weight. Also optionally available in a latex-free version (please enter LF for latex-free when ordering). High-quality precision sphygmomanometer-twin-tube model with black lacquered metal casing.

- › Life-long precision thanks to practically non-ageing and pressure-resistant measuring system up to 600 mm Hg.
- › Stainless metal clamp for hooking the pressure gage to the cuff.
- › Precision air release valve with fine adjustment and wear-free
- › Microfilter protects air release valve and measuring system.design for better readability of values.
- › Linear scale ø 49 mm from 0 to 300 mm Hg.
- › Latex bulb for pumping up the cuff.
- › Maximum error tolerance +/- 3 mm Hg.
- › In vinyl zipper case.
- › Comprehensive selection of different sizes and versions of cuffs.

No. 1350-100	exacta® manometer w/o cuff
No. 1350-102	exacta®, bandage cuff (adult)
No. 1350-105	exacta®, hook cuff (adult)

No. 1350	exacta®, velcro cuff (adult)
No. 1350-130	exacta®, velcro cuff (children)
No. 1350-123	exacta®, velcro cuff (large adult)

Please visit riester.de for more information