

VOLUMETRIC INCENTIVE SPIROMETER



- A Inspiratory indicator** encourages patients to maintain an appropriate rate of inspiration, when using the incentive spirometer the patient should aim to keep the yellow ball in the smiley face section.
- B Adjustable target indicator** provides a volumetric guideline.
- C Volume scale printed on both sides** dual-sided, clear identification of volumes.

The volumetric incentive spirometer helps patients perform and monitor their breathing exercises. The patients are encouraged to breathe in through the device as slowly and deeply as possible. The volumetric incentive spirometer provides visual guidance of the volume achieved during inspiration whilst encouraging the patient to maintain an appropriate inspiratory flow rate, by ensuring the yellow ball is kept within the smiley face section of the device.



- Adjustable target indicator provides a volumetric guideline for adequate lung expansion.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.



Order Information

Code	Description	Units Per Case
TB-93250	Volumetric Incentive Spirometer, 2500ml	48
TB-93500	Volumetric Incentive Spirometer, 5000ml	48



BESMED HEALTH BUSINESS CORP.

No. 2, Lane 106, Wu-Kong 3rd Rd, Wu-Ku District, New Taipei City, Taiwan. 24889

Tel: +886-2-2290-3959

Fax: +886-2-2299-9076

E-mail: info@besmed.com

Web-site: www.besmed.com



2021-02 Rev.0