



Dining Table

MWA-200

Table size: L765 × W385 mm
 Adjustable height: 720 ~ 1120 mm
 Frame: Powder coating / Electroplate with H shape base
 Mechanical system for height adjustable
 Table top: Wooden laminate with higher rim
 K/D design

Options:

- MWA-200-1 Wear-resist cover (MOQ request)
- MWA-200-2 Tilt table top
- MWA-200-3 2 sections table
- MWA-200-4 Drawer



MWA-200A

Table size: L900 × W400 mm
 Adjustable height: 625 ~ 925 mm
 Base: Plastic material
 Gas-spring system for height adjustable
 Table top: Wooden laminate with higher rim
 K/D design

Options:

- MWA-200A-1 Plastic cover
- MWA-200A-2 Higher table (MOQ request)



Recline Chair

MWA-300

Chair size: L800 × W600 × H950 mm
 Bed size: L1850 × W600 × H700 mm
 Frame: Stainless steel / powder coating steel
 Artificial leather with PU foam



Home Care Chair

MWA-400

Dimension: H1200 × W770 × D960 cm
 Adjustable back rest and foot rest
 Ideal for use in dialysis, blood collection, respiratory and cardiac care
 Side panels "Pop-off" for easy cleaning and maintenance
 5" Casters (2 with locks) make it easy to move and secure recliner
 Loading capacity: 250 lbs



Recline Chair

MWA-500

Chair size: L750 × W690 × H1050 mm
 Bed size: L1950 × W690 × H450 mm
 Frame: Powder coating steel (silver or white)
 Artificial leather with PU foam

Foldable arm rest
 Safety handle for extend function
 2" direction wheel × 4
 3" brake caster × 2



